

SPECIFICALLY DESIGNED FOR FEEDING AND DECOMPRESSION


Halyard Health, a quality leader in the enteral access market and the first to develop a gastrostomy tube specifically designed for feeding and decompression, offers the HALYARD* G Feeding Tube, designed to deliver optimal performance and value.

The unique design of the HALYARD* G Feeding Tube extends the balloon beyond the distal tip at recommended fill volumes, thereby decreasing the potential for gastric wall irritation or erosion. As well, the tip is tapered, easing catheter insertion and enhancing patient comfort, and a Luer-Lok® balloon port connector makes it easier for one-handed operation.



ADDITIONAL FEATURES OF THE HALYARD* G FEEDING TUBE INCLUDE:

- **High clarity medical grade silicone** construction for visibility and drapability
- **Clearly marked feeding, medication and balloon ports** for ease of identification
- **Ventilated external retention ring** designed to allow air to circulate around the stoma site and to reduce pressure 360° from the stoma
- **Radiopaque stripe** to aid in catheter visualization
- **Graduated centimeter markings** for tube position verification
- **Luer-Lok balloon port** to ensure secure syringe connection (also compatible with Luer-slip syringes)
- **Color-coded balloon port** according to catheter French size for ease of identification.




HALYARD* G FEEDING TUBE

HALYARD* G FEEDING TUBES



Stock #	Outer Diameter Size	Balloon Volume	Packaging
● 0112-14LV	14Fr	3-5 ml Balloon	1 Each
● 0112-16LV	16Fr	3-5 ml Balloon	1 Each
● 0112-18	18Fr	7-10 ml Balloon	1 Each
● 0112-20	20Fr	7-10 ml Balloon	1 Each
● 0112-22	22Fr	7-10 ml Balloon	1 Each
● 0112-24	24Fr	7-10 ml Balloon	1 Each



For more information, please visit:

www.halyardhealth.com/digestivehealth

Call 1-844-HALYARD (1-844-425-9273)
in the United States and Canada.

www.halyardhealth.com

