

Isolation Gown

With Thumb Hooks

ENCOURAGES COMPLIANCE AND ENHANCES PROTECTION

Our open-backed Over-the-Head Isolation Gown with Thumb Hooks features an easy on/easy off design that allows you to don and doff it in just seconds. It's generously sized for freedom of movement and made of cool, comfortable tri-layer fabric. And the open back provides ventilation for added comfort.

Intended for tasks involving light fluid exposure, this gown also includes thumb hooks for enhanced protection. When you slide your thumbs into the openings, your sleeves are held securely in place while you don your gloves. No sleeve ride-up. No glove slippage. Your hands and wrists remain covered until you're ready to remove your protective apparel.



Enhanced Features:




**Fast, easy-don/easy-doff
over-the-head design**



**Thumb hooks help keep
gloves in place while you
work**



**Open back provides
increased ventilation and
comfort**



The HALYARD* Over-the-Head Isolation Gown with Thumb Hooks is Made with Protective Tri-Layer SMS Fabric

Outer layers—Spunbond polypropylene material that is cloth-like, yet tough and abrasion-resistant

Middle layer—Intricate meltblown polypropylene micro fibers that filter out many water-based liquids

HALYARD* TRI-LAYER OVER-THE-HEAD ISOLATION GOWN WITH THUMB HOOKS

Code	Description	Size	Fabric Weight	Gowns/Case	Bags/Case
29549	Thumb Hooks, Yellow, Perforated Open Back	Universal	Medium	100 Case	10
29516	Thumb Hooks, Yellow, Perforated Open Back	XX-Large	Medium	100 Case	10
32574	Thumb Hooks, Yellow, Perforated Open Back, 11" x 5.5" Fold	Universal	Medium	120 Case	12
32572	Thumb Hooks, Yellow, Perforated Open Back, 11" x 5.5" Fold	XX-Large	Medium	120 Case	12

HALYARD* BASICS* TRI-LAYER OVER-THE-HEAD ISOLATION GOWN WITH THUMB HOOKS

Code	Description	Size	Fabric Weight	Gowns/Case	Bags/Case
43148	Thumb Hooks, Yellow, Perforated Open Back	Universal	Light	150 Case	6
43149	Thumb Hooks, Yellow, Perforated Open Back	X-Large	Light	150 Case	6

Gown is quick and easy to doff



KNOWLEDGE NETWORK* Accredited Education
In-Service Training and Technical Support
Credentialed Sales Representatives
Tools & Best Practices
Clinical Research
Commitment to Excellence

www.halyardhealth.com



*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2014 HYH. All rights reserved.

IP52 C14167