Sensitive Skin Comfort

MALYARD

Surgical and Procedure Masks Fluid-Resistant Protection

THE MASKS YOU'VE BEEN WAITING FOR ARE HERE

Halyard has always offered the best in facial and respiratory protection. Now, with our new line of Level 1 fluid-resistant face masks designed for sensitive skin, you can have the comfort you've always wanted, without sacrificing protection.

Introducing three new Level 1 Surgical and Procedure Masks.

- Feature an inner lining made of our exclusive SO SOFT* fabric, a unique bicomponent material that was rated softest and silkiest among leading masks tested.1
- Provide ASTM F2100-11 Level 1 protection, for procedures with risk of low fluid exposure.
- Dye-free inner lining, designed specifically for those with sensitive skin.



For assured Level 1, 2, or 3 fluid protection.

Look for the symbol on the package!

I FVFI 1	

Test Method	Test Description	ASTM F2100-11	Masks
ASTM F1862	Fluid Resistance with synthetic blood, in mm Hg	80 mm Hg	80 mm Hg
MIL-M-36954C	Delta P Differential pressure, mm H ₂ O/cm ²	<4.0 mm H₂O	2.7
ASTM F2101	Bacterial Filtration Efficiency (BFE), %	≥ 95%	99.9%
ASTM F2299	Submicron Particle Filtration (PFE) @ 0.1 μ, %	≥ 95%	98.8%
16 CFR Part 1610	Flammability	Class 1	Class 1

† Based on the calculated median.



Halyard Level 1 Surgical and Procedure Masks feature the silky feel of our dye-free SO SOFT* inner lining designed for sensitive skin. Offers ASTM F2100-11 Level 1 fluid protection, for low fluid exposure situations.





41803



HALYARD* SURGICAL AND PROCEDURE MASKS							
Code	Description	Masks/Box	Boxes/Case	Case Count			
41802	Fog-Free Procedure Mask, with SO SOFT* Lining, Earloops, White, ASTM Level 1 Fluid-Resistant	50	10	500			
41803	Fog-Free Procedure Mask, with WrapAround Visor, SO SOFT* Lining, Earloops, White, ASTM Level 1 Fluid-Resistant	25	4	100			
41805	Surgical Mask, with SO SOFT* Lining, Ties, White, ASTM Level 1 Fluid-Resistant	50	6	300			

TRY OUR COMPLETE LINE OF ASTM F2100-11 LEVEL 1 FLUID-RESISTANT MASKS							
Code	Description	Masks/Box	Boxes/Case	Case Count			
25867	Procedure Mask with SO SOFT* Lining, Earloops, Yellow	50	10	500			
25868	Procedure Mask with SO SOFT* Lining, Earloops, Lavender	50	10	500			
25869	Procedure Mask with SO SOFT* Lining, Earloops, Blue	50	10	500			
28802	Surgical Mask with SO SOFT* Lining, Ties, Blue	50	6	300			
28806	Fog-Free Surgical Mask with SO SOFT* Lining, Ties, Green	50	6	300			

For more information, visit our website at www.halyardhealth.com, or contact your Halyard representative.



KNOWLEDGE NETWORK* Accredited Education In-Service Training and Technical Support Credentialed Sales Representatives Tools & Best Practices Clinical Research Commitment to Excellence

