



# ORTHOPEDIC DRAPE

## DRAPE DESIGNED TO PROTECT YOU AND YOUR PATIENTS IN ORTHOPEDIC PROCEDURES


Now you have more choice than ever with HALYARD\* Orthopedic Drapes and Packs that deliver performance, protection and convenience to meet the draping needs of orthopedic surgery.

## THE RIGHT DRAPE FOR EVERY PROCEDURE

As a recognized leader in the development of high-performance medical fabrics, Halyard has rights to over 600 patents for advanced nonwovens materials.

HALYARD\* Surgical Drapes made with SMS and CONTROL\* Plus reinforcement fabric provide the critical balance of properties your demanding work requires.

- CONTROL\* Plus reinforcement fabric provides the highest rating for flame resistance<sup>1</sup> and resistance to laser ignition<sup>2</sup>.
- Low-lint<sup>3</sup> and abrasion-resistant fabric<sup>4</sup>.
- Fabric resistant to tearing<sup>5</sup> and strikethrough<sup>6</sup>.
- Surgical Drapes with CONTROL\* Plus reinforcement fabric meets AAMI Level 4 liquid barrier performance test requirements<sup>6</sup>.




# INTRODUCING OUR NEW ORTHOPEDIC DRAPE

## ORTHOARTS\* AND EXTRA LONG LOWER EXTREMITY DRAPE

Our extremity drapes provide a full array of features to meet the special needs of orthopedic surgery, including:

- Extra length for coverage where it's needed
- Conformable 2.5" diameter fenestration fits around patient's limb
- Absorbent CONTROL\* Plus reinforcement fabric
- Convenient, built-in armboard covers and tube holders


Our OrthoARTS\* Extremity Drapes are made with heavy-duty SMS fabric to withstand the rigors of orthopedic procedures.



Code	Description
<b>29219</b>	<b>OrthoARTS* Extra Long Lower Extremity Drape</b>
29220NS	OrthoARTS* Extra Long Lower Extremity Drape

All codes in bold indicate a sterile product  
Non-sterile product available for custom trays.

## EXTRA LONG LOWER EXTREMITY DRAPE



Code	Description
<b>29222</b>	<b>Extra Long Lower Extremity Drape</b>
29223NS	Extra Long Lower Extremity Drape

All codes in bold indicate a sterile product  
Non-sterile product available for custom trays.


## UNIVERSAL SPINE DRAPE

This latest addition to our orthopedic procedure line features:

- 4" x 12" fenestration, with generous clear instrument pouches on both sides
- Absorbent CONTROL\* Plus fabric reinforcement
- Convenient, built-in armboard covers and hook-n-loop tube holders

Code	Description
<b>29221</b>	<b>Universal Spine Drape</b>
29224NS	Universal Spine Drape

All codes in bold indicate a sterile product  
Non-sterile product available for custom trays.



## SHOULDER SPLIT PACK

Our Shoulder Pack includes three popular drapes available only in this pack.

- Shoulder split drape, 100" x 62", with adhesive circular fenestration
- Impervious split drape, 76" x 60" with a 6" split
- Body split drape with 10" x 24" split

## OTHER PACK COMPONENTS:

- XL and L ULTRA Surgical Gowns
- Mayo Stand Cover
- Back Table Cover, 50 in x 90 in
- 4 Utility Drapes
- Suture Bag

## Code Description

**29249** **Shoulder Split Pack**

All codes in bold indicate a sterile product  
Non-sterile product available for custom trays.

# ORTHOPEDIC DRAPES

5 Digit Code	Type	Description	Case Count
<b>29222</b>	<b>Sterile</b>	<b>Extra Long Lower Extremity Drape</b>	<b>13</b>
29223NS	Non-sterile	Extra Long Lower Extremity Drape	14
<b>29219</b>	<b>Sterile</b>	<b>OrthoARTS* Extra Long Lower Extremity Drape</b>	<b>8</b>
29220NS	Non-sterile	OrthoARTS* Extra Long Lower Extremity Drape	8
<b>29221</b>	<b>Sterile</b>	<b>Universal Spine Drape</b>	<b>8</b>
29224NS	Non-sterile	Universal Spine Drape	9
<b>29249</b>	<b>Sterile</b>	<b>Shoulder Split Pack</b>	<b>4</b>

All codes in bold indicate a sterile product. Non-sterile product available for custom trays.



KNOWLEDGE NETWORK\* Accredited Education  
In-Service Training and Technical Support  
Credentialed Sales representatives  
Tools & Best Practices  
Clinical Research  
Commitment to Excellence

For more information, please visit:

[www.halyardhealth.com](http://www.halyardhealth.com)

Call 1-844-HALYARD (1-844-425-9273)  
in the United States and Canada.



\*Registered Trademark or Trademark of Halyard Health, Inc.  
© 2015 HYH. All rights reserved. C15547

## Ask about our C-Arm and Banded Bag!

- 1 Flammability Standard: 16FR 1610: 2003-Class 1
- 2 ISO 11810: 2005 Part One
- 3 Test method ISO 9073-10; Criteria <4.0 Log10 (lint count).
- 4 Test method ASTM D4970:2007 Criteria  $\geq 3$
- 5 Test method DSS87-08:2008
- 6 Test method ASTM F1670: 2008